



TECHNICAL DATA

MASKS



info@virostop.be | www.virostop.be | +32 4 330 18 60



TYPE 1 WASHABLE FOLD BARRIER MASKS

Barrier Masks Non Sanitary Use (Cat 1) for people with occasional contact with the public, not suitable for nursing staff in contact with patients.

Our range of textile masks produced in our factories meets European standards.

100% Cotton made, usable 8 hours and washable 50 times at 60°C .

Our elastics, according to OEKO-TEX standards, support a very large number of washings.



We offer ear or head masks, in 4 sizes (3 child sizes and 1 adult size), masks with interchangeable filter, or even waterproof (see details on www.virostop.be).

VIROSTOP masks are customizable.

TECHNICAL DATA

Bacterial filtration efficiency (BFE)	_____	≥95%
Differential pressure (Pa/cm ²)	_____	125
Spray resistance pressure	_____	70
Microbial cleanliness (cfu/g)	_____	not measured by DGA
Use	_____	8 hours
Washable	_____	50x
Main material	_____	100 % cotton

ECO-FRIENDLY

The VIROSTOP mask is an intelligent solution, its design makes it effective, comfortable, durable 50 x 8 hours (which allows among other things to amortize your investment/ advertising when printed) but also eco-friendly.

After its service, our mask made of 100% cotton, will end its life in a textile recycling industry

This choice has **2 main advantages**:

- 1 it respects the **balance of the planet**
- 2 it testifies to your **involvement** in the **protection of the environment**

ADVICE

- Store in a clean, dry and dark place
- This mask is strictly adapted to the pandemic COVID 19.
- Social distancing must be applied despite wearing the mask.



Rue de l'Environnement 22
4100 SERAING - BELGIQUE
Tel: +32 4 330 18 60
E-mail: info@contipharma.com

AND ALSO :
OUR
GLOVES

ALCOHOL-FREE NATURAL
DESINFECTANT
GEL

OUR
GEL

