



**CONTIPHARMA**  
a member of Bedelco group

**COVID19**

## Disposable Civil and Medical masks



### FEATURE

(1) Disposable mask blue colour, with earloop. 3 layer of nonwoven fabric.

1st. Layer: 25 gm. spunbond nonwovens

2nd. Layer: 25 gsm Meltbond nonwovens

3rd. Layer: 22 gsm. spunbond nonwovens

Package: 50pcs/box, 2000pcs/carton

Daily capacity : 50,000 PCS/day

(2) Disposable mask white colour, with earloop, 3 layer of nonwoven fabric

1st. Layer: 25 gm. spunbond nonwovens

2nd. Layer: 25 gsm Meltbond nonwovens

3rd. Layer: 22 gsm. spunbond nonwovens

Package: 50pcs/box, 2000pcs/carton

Daily capacity : 1000,000 PCS/day

---

### CERTIFICATE

Examination Intent :

Examination the completeness of the Technical Documentation according to the requirements of Annex II excluding section 4 (Module H) of 93/42/EEC

Harmonised standards :

EN 13060:2014; EN 60601-1: 2006; EN ISO 13485: 2016; EN ISO 15223-1:2016

Directive 2016/425 PPE



# CONTIPHARMA

a member of Bedelco group

## COVID19

### DISPOSABLE FACE MASK - BKDSFM



### INTRODUCTION

It is widely used in ordinary environment; to block the mouth and nasal cavity from exhaling or spraying pollutants by covering the user's mouth, nose and jaw.

### FEATURES

- Moderate tightness provide coziness for ears.
- Light and very low resistance to breathing.
- Three-layer anti-wrinkle design.

### TECHNICAL PARAMETERS

|                                 |   |
|---------------------------------|---|
| Model                           | BK-DSFM   |
| Type                            | Flat earloop type   |
| Structure                       | Consisting of mask body, mask band and nose clip  |
| Material                        | PP non-woven fabric, rubber band  |
| Instruction                     | Turn folded side out, nose bar side up  |
| Three layers material           | The modeling layer is spundonted non-woven fabric<br>The filter layer is melt-blown non-woven fabric<br>The affinity layer is spundonted non-woven fabric |
| Product protection effect level | Level A   |
| External Size                   | 180*95mm (with 3 layers)  |
| Recommended using time          | 4 ~ 6 hours   |
| Validity                        | 2 years   |
| Storage condition               | Keep in a sealed, cool and dry condition  |
| Packing size                    | 420*300*590mm (50 pcs/bag, 2500pcs/Ctn)   |
| Gross weight                    | 7.5 kg  |



# CONTIPHARMA

a member of Bedelco group

## COVID19

### DISPOSABLE CIVIL FACE MASK



Product

Package : Quantity : 50pcs/BOX - 2000pcs/carton  
Carton size: 530\*400\*400mm,  
Weight : 9 kg/Ctn.  
Standard: GB2626-2006  
Certification : CE Approved  
MOQ: 100,000 PCS  
Lead time: 3-5 days  
Less or Large numbers need double check.  
1 million 12-15 days

### DISPOSABLE MEDICAL FACE MASK



Product

Package : Quantity : 50pcs/BOX - 2000pcs/carton  
Carton size: 600\*440\*425mm,  
Weight : 10,45 kg/Ctn.  
Standard: GB2626-2006  
Certification : CE and FDA Approved  
MOQ: 100,000 PCS  
Lead time: 3-5 days  
Less or Large numbers need double check.  
1 million 12-15 days





# CONTIPHARMA

a member of Bedelco group

## COVID19

### DISPOSABLE MEDICAL SURGICAL FACE MASK



Product

Package : Quantity : 50pcs/BOX - 2000pcs/carton  
Carton size: 580\*380\*300mm,  
Weight : 9 kg/Ctn.  
Standard: GB2626-2006  
Certification : CE Approved  
MOQ: 100,000 PCS  
Lead time: 3-5 days  
Less or Large numbers need double check.  
1 million 12-15 days

### KN95 FFP2 FACE MASK



Product

Package : Quantity : 1pc/plastic bag - 10 pcs/box - 2000pcs/Ctn  
Carton size: 700\*625\*540mm,  
Weight : 18 kg/Ctn.  
Standard: GB2626-2006  
Certification : CE Approved  
MOQ: 100,000 PCS  
Lead time: 12 days  
Less or Large numbers need double check.

#### KN95 FFP2 FACE MASK



Product

Package : Quantity : 2pcs/Plastic bag - 400bags/Ctn - 800 pcs/ctn  
Carton size: 580\*400\*520mm,  
Weight : 9,8 kg/Ctn.  
Standard: GB2626-2006  
Certication : CE Standard : EN149:2001+A1:2009  
CE and FDA Approved  
MOQ: 100,000 PCS  
Lead time: 10 days for 500,000 pcs  
Less or Large numbers need double check.

#### KN95 FFP2 CHILDREN FACE MASK



Product

Package : Quantity : 10 pcs/box - 420 pcs/Ctn  
Carton size: 490\*330\*395mm,  
Weight : 4,3 kg/Ctn.  
Standard: GB2626-2006  
Certification : CE and FDA Approved  
Lead time: 7 days for 200,000 pcs  
Less or Large numbers need double check.



# CONTIPHARMA

a member of Bedelco group

## COVID19

### N95 Face mask



#### Product

Package : Quantity : 20pcs/box - 24 box/Ctn - 480 pcs/Ctn  
Carton size: 620\*290\*290mm,  
Weight : 4,3 kg/Ctn.  
Certification : CE and FDA approved  
MOQ: 10,000 Pcs  
12 days for 100,000 pcs  
20 days for 200,000 pcs  
Less or Large numbers need double check.

CONTIPHARMA

Rue de l'Environnement 22  
4100 Seraing  
Belgium

Mobile : +32 475 23 54 21

Email : [info@contipharma.com](mailto:info@contipharma.com)

Web : [www.contipharma.com](http://www.contipharma.com)